

# Thermo Scientific Dionex Corona Veo Charged Aerosol Detector

Discover What You're Missing

We are dedicated to providing innovative solutions for analytical laboratories. To meet these demands we deliver premier instrumentation across our complete portfolio, including new products designed to deliver top performance and productivity for today's most difficult chromatography challenges. Our instruments and automation systems provide the answers that allow you to advance your work and achieve success.



Corona Veo Charged Aerosol Detector.



UltiMate 3000 LC system with Corona Veo detector.

The Thermo Scientific™ Dionex™ Corona™ Veo™ detector represents the evolutionary refinements in instrumentation design that uses the widely acknowledged charged aerosol detection (CAD) technology for liquid chromatography; including traditional HPLC, UHPLC, and micro LC.

The Corona Veo charged aerosol detector delivers performance that other LC detectors simply cannot match. Consistent response independent of analyte chemical structure is a unique characteristic of charged aerosol detection, which lets you estimate relative amounts even without standards. Charged aerosol detection technology helps you see analytes that other systems fail to detect. Any nonvolatile and most semi volatile analytes with or without a chromophore can be measured using this technology. With a predictable response, the Corona Veo detector will quickly become your first choice for HPLC and UHPLC analyses.

## Features

- Near universal detection
- Fast data sampling ensures full compatibility with UHPLC
- Improved semi-volatile sensitivity with the new FocusJet™ concentric nebulizer
- Adjustable evaporation temperature for improved analyte response
- Intuitive operation
- Integration with any HPLC/UHPLC system
- Available software drivers for Thermo Scientific™ Dionex™ Chromeleon™, Empower® 2 and 3, ChemStation®, EZChrom® and OpenLab® Chromatography Data System software.

**UHPLC<sup>+</sup>**  
focused

**Thermo**  
SCIENTIFIC

The Corona Veo detector can be used with the most up-to-date UHPLC technology, such as Thermo Scientific™ Dionex™ UltiMate™ 3000 LC systems, combining all the benefits of charged aerosol detection with the high speed and increased resolution of UHPLC. Whichever configuration you choose, you will get a highly integrated solution with optimized fluidic connections and single-point intelligent control through Chromeleon CDS software.

<b>SPECIFICATIONS</b>		
<b>Description</b>	<b>Corona Veo</b>	<b>Corona Veo RS</b>
Operating Mode	Charged Aerosol Detection	
Nebulization	FocusJet™ Concentric flow design	
Mobile Phase Flow Rate	0.2–2.0 mL/min	0.01–2.0 mL/min
Wetted Materials	stainless steel (type 316), Nitronic® 60 stainless steel, PEEK™, Simriz®, aluminum, fused silica, and PTFE: Valcon H, Valcon E (Corona Veo RS only)	
Digital Data Collection Rate (max.)	100 Hz	200 Hz
Digital Noise Filtering	4th order low-pass Bessel	
Optional Analog Signal Output	0–1 V DC (field installed)	
Full Scale Analog Output Range	1 pA to 500 pA in 1-2-5 sequence	
Filter Time Constants	Selectable in numerical sequence (1-2-5 increments)	
Standalone Control Interface	Integrated color LCD touch screen	
Evaporation Temperature	Select: 35 °C or 50 °C	Settable Range: Ambient + 5 °C to +100 °C
Warm-Up Time	< 30 min to 35 °C evaporation T	
Integrated Stream Switching	Not available	TTL controlled, 6-port, 2-position valve
Inlet Gas Supply	Compressed Air or Nitrogen Inlet Pressure 70–80 psig (482–551 kPa)	
Gas Pressure Control	Manual	Electronic
PC Connection	Available detector functions controllable via USB 2.0 through Chromeleon data system (cable provided)	
Rear Interface	AC plug, power switch, USB port, Digital I/O (four TTL Inputs, two Relay Outputs)	
Dimensions (h × w × d)	22.9 cm × 44.5 cm × 55.9 cm (9 in. × 17.5 in. × 22 in.)	
Weight	14.3 kg (31.5 lbs.)	
Power Requirements	100/240 VAC, 50/60 Hz, 100 VA	
Environmental Temperature Range	15 °C to 35 °C at 12% to 80% RH, non-condensing	
Safety Certifications	UL/IEC 61010-1, 3 <sup>rd</sup> Edition FCC Part 15/ICES-003 class A	

## Ordering Information

Detector Modules	Part Number
Corona Veo Charged Aerosol Detector	5081.0010
Corona Veo RS Charged Aerosol Detector with electronic gas regulation and automated stream-switching	5081.0020

  

Options	Part Number
Analog signal output kit for Corona Veo detector	6081.0010
Adjustable analytical flow splitter (1:1 to 1:20)	70-6337

### [www.thermoscientific.com/Veo](http://www.thermoscientific.com/Veo)

©2013 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. ChemStation, EZChrom and OpenLab are registered trademarks of Agilent Technologies. Empower is a registered trademark of Waters Corporation. PEEK is a registered trademark of Victrex PLC. Nitronic is a registered trademark of AK Steel Corp. Simriz is a registered trademark of Freudenberg-NOK. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales



Thermo Fisher Scientific, Sunnyvale, CA  
USA is ISO 9001:2008 Certified.

<b>Australia</b> +61 3 9757 4486	<b>Denmark</b> +45 70 23 62 60	<b>Japan</b> +81 6 6885 1213	<b>Switzerland</b> +41 62 205 9966
<b>Austria</b> +43 1 333 50 34 0	<b>France</b> +33 1 60 92 48 00	<b>Korea</b> +82 2 3420 8600	<b>Taiwan</b> +886 2 8751 6655
<b>Belgium</b> +32 53 73 42 41	<b>Germany</b> +49 6126 991 0	<b>Netherlands</b> +31 76 579 55 55	<b>UK/Ireland</b> +44 1442 233555
<b>Brazil</b> +55 11 3731 5140	<b>India</b> +91 22 6742 9494	<b>Singapore</b> +65 6289 1190	<b>USA and Canada</b> +847 295 7500
<b>China</b> +852 2428 3282	<b>Italy</b> +39 02 51 62 1267	<b>Sweden</b> +46 8 473 3380	

# Thermo

SCIENTIFIC

Part of Thermo Fisher Scientific