

建宏層析企業股份有限公司

VERCOTECH INC.,

地址：台北市內湖區港墘路 221 巷 19 號 5 樓

電話：02-2657-8289 傳真：02-2657-8389 Email: vercopak@vercopak.com.tw






Cadenza HS-C18




World's First Reproducible Direct Injection Column

- 不需萃取即可直接注入分析，節省時間
- 立刻溶洗出分子量較大的 Protein (>30kDa)，接著再分離分子量較小的化合物
- 革命性的技術，達成不同批次絕佳的再現性
- 最適合分析疏水性化合物

Hybrid ODS 固定相中，分子量與 Protein 的排除間的關係



● Bovine Serum 的直接注入分析



● Ordering Information for Cadenza HS-C18 Columns

3µm Column, Pressure limits of up to: 20MPa, 250 bar, 3,000 psi

Column Length	Internal Diameter						
	1.0 mm	1.5 mm	2.0 mm	3.0 mm	4.6 mm	6.0 mm	10.0 mm
10			HS020	HS030	HS000		
20			HS029	HS039	HS009		
30	HS011	HS071	HS021	HS031	HS001	HS061	HS0P1
50	HS012	HS072	HS022	HS032	HS002	HS062	HS0P2
75	HS013	HS073	HS023	HS033	HS003	HS063	HS0P3
100	HS014	HS074	HS024	HS034	HS004	HS064	HS0P4
150	HS015	HS075	HS025	HS035	HS005	HS065	HS0P5
250	HS016	HS076	HS026	HS036	HS006	HS066	HS0P6
500					HS007		

Guard Column System for Cadenza HS-C18

	Internal Diameter						
	1.0 mm	1.5 mm	2.0 mm	3.0 mm	4.6 mm	6.0 mm	10.0 mm
Guard Holder	GCH01S	GCH01S	GCH01S	GCH01S	GCH01S	GCH01S	GCH02M
Guard Cartridge (Set of 3)	GCHS0C	GCHS0C	GCHS0S	GCHS0S	GCHS0S	GCHS0S	GCHS0M

All of our stationary phases can also be made in the following internal diameters:
Nano: 0.05mm, 0.075mm **Capillary:** 0.1mm, 0.3mm, 0.5mm